

LabVIEW Guide – Installation Instructions

LabVIEW source available from: https://github.com/D-TACQ/acq400_lv/releases/tag/2.0.0

The screenshot shows the GitHub repository page for `D-TACQ / acq400_lv`. The repository has 1 Watch, 0 Stars, and 2 Forks. The `Releases` tab is selected, showing the latest release `Full_House_includes_ACQ1014` by `petermilne`, released 20 hours ago. The release is verified and has a commit hash of `038c2ee`. Under the `Assets` section, there are two download links: `Source code (zip)` and `Source code (tar.gz)`. A yellow callout bubble points to the `Source code (zip)` link with the text: "Download source code as a zip file and extract to a known location. All paths given in this document are relative to the install location of the LabVIEW source code. Any release from 2.0.0 onwards will be suitable."

Search or jump to... [Pull requests](#) [Issues](#) [Marketplace](#) [Explore](#)

D-TACQ / `acq400_lv` [Watch](#) 1 [Star](#) 0 [Fork](#) 2

[Code](#) [Issues](#) 0 [Pull requests](#) 0 [Projects](#) 0 [Wiki](#) [Insights](#)

Releases [Tags](#)

Latest release

2.0.0
038c2ee
Verified

Full_House_includes_ACQ1014

petermilne released this 20 hours ago

Assets

- [Source code \(zip\)](#)
- [Source code \(tar.gz\)](#)


Merge pull request #4 from seanalsop/master


Working transient capture and acq1014 plotting tools.

Download source code as a zip file and extract to a known location. All paths given in this document are relative to the install location of the LabVIEW source code. Any release from 2.0.0 onwards will be suitable.

LabVIEW Guide – Transient capture

Vi available from: acq400_1\Source Code\Combined Status and Transient\Combine Status Monitor and Transient Capture (Single UUT).vi







After choosing a configuration file use the soft trigger button to configure a soft trigger unless you want to use an external trigger with your UUT.

Use the configuration file path button to view the list of .ini files. Choose a single file and select OK.

Press the “run vi” button in the top left of the screen and the graph and status monitor will be ready to display data as shown above.



Sub Panel



- ☒ CH01
- ☒ CH02
- ☒ CH03
- ☐ CH04
- ☐ CH05
- ☐ CH06
- ☐ CH07
- ☐ CH08



Press the "FIRE" button and the UUT is armed.

If the soft trigger button on the right is selected the UUT will soft trigger and the data will be displayed in the graph above.

Sub Panel

Device State
POST_PROCESS ☒ Demuxing

Pre-Buffered Samples
0

Post-Buffered Samples
100000

Total Samples
262144

Configuration file path
[Browse]

Stop

configuration file path
[Browse]

Post-Trigger Samples
100000

Soft Trigger
☒

UUT Name
ACQ1001_329


stop
STOP

When the UUT is armed it will display its status here.

LabVIEW Guide – Hardware in the Loop

Vi available from: acq400_iv\Source Code\Hardware in the Loop (Advanced Demo)\Demo - Hardware in the Loop Application.vi

How the vi looks before it is run. Click the INI file path button and choose a single file for the HIL test.



Demo - Hardware in the Loop Application.vi Front Panel

File Edit View Project Operate Tools Window Help

15pt Application Font

Search

Open

<< acq400_lv > Source Code > Communications > Configuration

Search Configuration

Organize

New folder

Quick access

Desktop

Downloads

Documents

Pictures

Workspace

Lightshot

Combined Status ar

Configuration

LabVIEW

transient_with_statu

Dropbox

OneDrive

This PC

Network

Name

Date modified

Type

Size

ACQ1001_047.ini
ACQ1001_084.ini
ACQ1001_105.ini
ACQ1001_105_2.ini
ACQ1001_148.ini
ACQ1001_149.ini
ACQ1001_278.ini
ACQ1001_279.ini
ACQ1001_301.ini
ACQ1001_329 - Copy.ini
ACQ1001_329.ini
ACQ1001_329_example.ini
ACQ1001_336.ini
ACQ2106_054.ini
ACQ2106_081.ini
ACQ2106_087.ini
Multi UUT (Sean).ini
Multi UUTs (1014).ini

29/05/2018 14:19
29/05/2018 14:19
03/07/2018 11:29
03/07/2018 11:30
05/07/2018 16:09
05/07/2018 16:09
05/07/2018 12:33
05/07/2018 12:37
29/05/2018 14:19
12/07/2018 08:37
08/06/2018 14:03
12/07/2018 08:44
06/07/2018 10:37
29/05/2018 14:19
29/05/2018 14:19
08/06/2018 14:03
05/07/2018 15:12
05/07/2018 16:24

Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...
Configuration sett...

1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB
1 KB

File name:

Configuration File (*.ini)

OK

Cancel

- ☒ CH01
- ☒ CH02
- ☒ CH03
- ☒ CH04
- ☒ CH05
- ☒ CH06
- ☒ CH07
- ☒ CH08
- ☒ CH09
- ☒ CH10
- ☒ CH11
- ☒ CH12
- ☒ CH13
- ☒ CH14
- ☒ CH15
- ☒ CH16

Speed

0.5

Dynamic Channel

0

error ou


status

source

stop

STOP

If the default settings are okay for the hardware you are using then click the run button in the top left of the window and the loop test will begin.



LabVIEW Guide – Streaming data

Vi available from: acq400_iv\Source Code\Continuous Capture (Advanced Demo)\Demo - Continuous Capture.vi

Demo - Continuous Capture.vi Front Panel

File Edit View Project Operate Tools Window Help

15pt Application Font

Click the configuration file path button and choose a single ini file.

configuration file path

Buffer Block Size (Samples per channel)

256

Data Size

16 Bit (Signed Short)

Show Console

Waveform Graph

Amplitude

Time

Capture Rate (Bytes/Second)

0

stop


STOP


error out

status code

source

CH01 CH02 CH03 CH04 CH05 CH06 CH07 CH08 CH09 CH10 CH11 CH12 CH13





configuration file path

C:\Users\sean\

Buffer Block Size (Samples per channel)


16384

Data Size

16 Bit (Signed Short)

Show Console

Waveform Graph



To reduce data rate on the wire,
first reduce the number of active
channels using the bank mask
UUT command:

“set.site 1 bank_mask ABC”

which selects 12 channels, for example,
then modify the channel count
in the .ini file to suit

ACQ1001_329_example.ini - Notepad

File Edit Format View Help

```
[Carrier]
Name="acq1001_329"
IP=10.12.197.130
Port=4220
Timeout=2500
Type="Carrier"
DownloadContinuousDataPortMux=4210
DownloadCapturedDataPortMux=53000
ModuleSlots=2
```


```
[CarrierModule1]
Name="ACQ425ELF"
Port=4221
Timeout=2500
Type="AnalogIn"
Channels=12
```

```
[CarrierModule2]
Name="AO424ELF"
Port=4222
Timeout=2500
Type="AnalogOut"
Channels=16
```


LabVIEW Guide – acq1014 Interface

Vi available from: acq400_iv\Source Code\Combined Status and Transient\Run - Array of Status Monitors With Transient Captures.vi


Choose a configuration file using the configuration file path button.




Decide whether you wish to use a soft trigger or an external trigger to trigger your UUTs.









Once the correct ini file has been chosen run the vi from the button in the top left of the window.






Currently Selected Unit Under Test



The vi displays the status of all of the UUTs in the ini file. Use the buttons shown here to use the UUTs.

Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_278	IDLE		0	0	0
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_279	IDLE		0	200000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_148	IDLE		0	200000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_149	IDLE		0	200000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_278	IDLE		0	200000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_279	IDLE		0	200000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_148	IDLE		0	200000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_149	IDLE		0	200000	1310720



Currently Selected Unit Under Test

ACQ1001_278

Arm All

Graph All

Previous


Next

Arm One UUT

Graph One UUT

Here a single UUT has been graphed using the "Graph One UUT" button. Note that the graph window must be open before the UUT is armed and triggered in order to view the data from that capture.

Demo - Transient Capture (Multiple Shots)_for_sub_panel.vi:1330010 (clone)



The previous and next buttons can be used to choose which UUT is currently selected for display.

- ☒ CH01
- ☒ CH02
- ☒ CH03
- ☐ CH04
- ☐ CH05
- ☐ CH06
- ☐ CH07
- ☐ CH08

FIRE

Device Name

ACQ1001_278



[UUT]
MAIN
MENU

Currently Selected Unit Under Test

ACQ1001_278



Previous



Next

Arm All

Graph All

Abort

Arm One UUT

Graph One UUT

Close Graph

Device states are shown here when the UUTs are armed.

Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_278	POST_PROCESS		0	100000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_279	CLEANUP		0	100000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_148	POST_PROCESS		0	100000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_149	IDLE		0	100000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_278	POST_PROCESS		0	100000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_279	CLEANUP		0	100000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_148	POST_PROCESS		0	100000	1310720
Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_149	IDLE		0	100000	1310720

Currently Selected Unit Under Test

ACQ1001_278

Previous

Next

Arm All

Graph All

Abort

Arm One UUT

Graph One UUT

Close Graphs

Stop All

Device Name	Device State	Demuxing	Pre-Buffered Samples	Post-Buffered Samples	Total Samples
ACQ1001_278	IDLE		0		
ACQ1001_279	IDLE		0		
ACQ1001_148	IDLE		0		
ACQ1001_149	IDLE		0	100000	1310720
ACQ1001_278	IDLE		0	100000	1310720
ACQ1001_279	IDLE		0	100000	1310720
ACQ1001_148	IDLE		0	100000	1310720
ACQ1001_149	IDLE		0	100000	1310720

Once all tasks have been completed with the Interface use the "Stop All" button to stop all of the vi's from running.

LabVIEW Guide – ini files

Ini files available from: acq400_iv\Source Code\Communications\Configuration

C:\Users\sean\Documents\Workspace\LabVIEW\Sean\acq400_iv\Source Code\Communications\Configuration\ACQ1001_105.ini - Notepad++

File Edit Search View Encoding Language Settings Tools Macro Run Plugins Window ?

ACQ1001_105.ini

```
1 [Carrier]
2 Name="ACQ1001_105"
3 IP="10.12.197.26"
4 Port=4220
5 Timeout=2500
6 Type="Carrier"
7 DownloadContinuousDataPortMux=4210
8 DownloadCapturedDataPortMux=53000
9 ModuleSlots=1
10
11 [CarrierModule1]
12 Name="ACQ480ELF"
13 Port=4221
14 Timeout=2500
15 Type="AnalogIn"
16 Channels=8
17
```

MS ini file length : 269 lines : 17 Ln : 17 Col : 1 Sel : 0 | 0 Windows (CR LF) UTF-8 INS

Here is the structure of an ini file for use with the LabVIEW code. Adjust the values to reflect the correct settings for your UUT. IP address must be used instead of DNS name. Unless otherwise specified the ports should always be as stated here.

For every module that the UUT has it should have another section listed as such, using the CarrierModuleX nomenclature.